GENERAL CHEMISTRY LAB COMPONENT CHE101L FALL 2020 CONTENT: LAB 1



Dissolution Reactions: Heats of Dissociation

NAME	Shafayet	Jamil	Zim

STUDENT ID 1620346042

SECTION 01

NAME OF THE INSTRUCTOR

Md. Mainul Hossain

SIGNATURE & DATE

10/110/20

REPORT SUBMISSION DATE

EXPERIMENT 1

DISSOLUTION REACTIONS: HEATS OF DISSOCIATION

Heats (exothermic or endothermic) are associated with chemical reactions. Quantity of heat evolved or absorbed is directly proportional to the amount reacted. Consider the reaction:

Heat could be generated or absorbed in this reaction. When heat is generated/released from a chemical reaction it is called exothermic reaction (you can feel it by touching the reaction container (warmer) and when heat is absorbed the reaction is called endothermic (colder). When reactions occur in a reaction vessel (e.g., Beaker) in aqueous condition, formation and dissociation of chemical bonds occur simultaneously. Bond formation and dissociation involves heat energy of the system which is expressed by the term Q which is called enthalpy.

PROBLEM STATEMENT: Is heat energy related to chemical reactions, how?

This experiment is subdivided into two parts:

1. QUALITATIVE & II. QUANTITATIVE

PART I. QUALITATIVE

DATA COLLECTION

Place about 30 mL of distilled water into a 50 mL beaker. Suspend a thermometer (having 0.1°C division mark) into the beaker using thermometer clamp and ring stand. Please make sure that the thermometer is not touching the bottom of the beaker, as any movement of the beaker could break the thermometer. Record the temperature of water in the beaker in every 30 seconds for 240 seconds.




FIGURE 1: Experimental setup for dissolution reaction

Place a moderate amount (which would be 1 to 3 cm³) of supplied anhydrous magnesium sulfate (MgSO₄) to the beaker. Mix vigorously with the glass rod for 5 minutes. Record your observations. (2 points)

Mg100 -> 25°e-27°e -> Exothermic

Exothermie

Repeat this procedure with each of the following compounds: (2 points)

- a. Sodium Nitrate, NaNO3 > 250 24.50 > Endolherme
- b. Sodium Chloride, NaCl -> 24.75°c 24.5°c -> Endother mic
- c. Hydrated Calcium Chloride, CaCl2.2H2O -> 24.5°c- 25.25° > Exothermic
- d. Ammonium Nitrate, NH4NO3 -> 24.25°c 24°e -> Endotherme

DATA ANALYSIS:

what are the similarities and differences in the behavior of these compounds? Can you find out any generalization concerning all chemical reactions here? What conclusion can be drawn from these data? (4 points)

theats are associated with the chemical neartion. At first all compounds dissolved and chemical bond breaks. Energy is required to break the bonds. Then chemical reaction occurs and new bonds forms. Energy transforms and for some heat realeased and ton some heat absorbed. If total heat of energy of the system is greater than zero then the neartion is endothermic. If total heat of energy of the system the system is less than zero then the near the near than is endothermic.

PART II. QUANTITATIVE

DATA COLLECTION:

 Accurately weigh a 3 to 5 gm sample of MgSO₄ on the analytical balance. Record the exact mass here. For 4 different trials below measure four different weight samples (e.g., 1,2, 4 & 5 grams respectively).

b. Suspend the thermometer into a polystyrene cup/coffee cup. Make sure of the thermometer is not touching the bottom of the cup. Measure 20 mL of distilled water by a volumetric cylinder into the cup and stir for 240 second. Record the temperature in every 20 seconds. After 240 seconds add MgSO₄ with vigorous mixing while continuing to record data for 5 minutes.

c. Determine the temperature change, ΔT , for the reaction. This can be done from the difference of the highest temperature minus the slope of the line go through the points from first 240 seconds of data.

d. Draw a temperature vs. time graph. Draw the best curve through the points and point out what is happing in each part of the curve.

TRAILS

(I) Mass of MgSO ₄ 0.5 gm		(II) Mass of N	(II) Mass of MgSO ₄ 1.0 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)	
20	26.0	20	27.1	
40	27.3	40	28.0	
60	28.0	60	28.3	
80	30.0	80	29.3	
100	30.1	100	30.5	
120	30.5	120	31.5	
140	30.5	140	31.5	
160	30.5	160	31.5	
180	30.5	180	32.0	
200	29.8	200	31.8	
220	29.5	220	31.5	
240	29.1	240	31.0	
260	29.2	260	31.0	
280	28.8	280	30.8	
300	28.1	300	30.8	

(II <u>I</u>) Mass of MgSO ₄ 1 <u>.5 gm</u>		(IV) Mass of MgSO ₄ 2.0 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)
20	25.5	20	27.5
40	25 - 5	40	28.5
60	25.5	60	29.5
80	27	80	30.5
100	29	100	31.5
120	32	120	32.5
140	33.25	140	32.5
160	33.25	160	32.5
180	33	180	32.3
200	33	200	31.5
220	33	220	31.5
240	32.75	240	31.0
260	32.5	260	31.0
280	22.25	200	

Nater

DATA ANALYSIS

- 1. What do you understand form the data you recorded and from the other trials? (4 points)

 avantity of the heat evolved on absorbed is

 dineetly proportional to the amount nearled.

 Heat changes with the time and goes to its

 wighest point. Then stanted does asing.
- Calculate the heat, Q & moles, n, of the reaction. Take help from the equitation Q = C x M x ΔT.
 Assume C = 4.18 Joules/gram °C and M is the mass of water (take the water density as 1.00 grams/cm³). (4 points)

Mgsou = 120 g

120g Mgsou
$$\rightarrow$$
 1 MOL


1.5 g Mgsou \rightarrow $\frac{1\times1}{120}$ MOL


1.5 g Mgsou \rightarrow $\frac{1'5}{120}$ MOL

= 0.0125 MOL

for 0.5 g mgsoy

Q = 376 5





MENTAL MODEL: Use the chemical equation given above to represent the dissolution reaction in this experiment. Draw a picture(s) which describes what is happening in atomic or in molecular level. How heat release or absorbed can be described from these pictures? (5 points)

total + (-) Keal/mol